



Model 6 Crimper

New OLED Screen Displays:

- Tool Size
- Battery Status
- Crimp Setting
- Cycle Result

Settings Menu Provides:

- Lifetime Statistics
- Recent Cycle Log
- 8 Language Options
- Reset Functions



Brushless Motor for Longer Life!



Electronic Crimping Tools



- Powerful brushless motor with a fast cycle and longer life
- OLED screen displays tool size, battery status, crimp setting and cycle result
- Settings menu shows lifetime statistics, recent cycle log, 8 language options and reset functions
- Long lasting battery with Lithium Ion technology
- Can be operated while charging
- Convenient controls on top of unit to adjust crimp
- For flip-off caps, available in 13 mm and 20 mm sizes*
- *Not for use with Flip-up caps*
- Optional base and mounting kit

CS electronic crimpers and decappers fit standard laboratory vials and seals with aluminum or two-part aluminum/steel caps. The electronically adjustable cycle is accurate, reproducible and can be adjusted with the controls on top of the unit. The tools are ergonomically designed to reduce strain and arm injury associated with many manual crimping tools. Available in 8 mm (crimper only), 11 mm, 13 mm, and 20 mm cap sizes.

Rechargeable: The tools can be recharged in one to two hours with the universal power supply included. The power supply includes a plug set for operation in most countries.

Description	Item No.
8 mm Electronic Crimper	029108-V6
11 mm Electronic Crimper	029111-V6
13 mm Electronic Crimper	029113-V6
20 mm Electronic Crimper	029120-V6
13 mm Electronic Crimper, for Flip Off Caps	029113-F0-V6
20 mm Electronic Crimper, for Flip Off Caps	029120-F0-V6

Description	Item No.
11 mm Electronic Decapper	029122-V6
13 mm Electronic Decapper	029123-V6
20 mm Electronic Decapper	029124-V6
Base with Mounting Kit	029158