



PRECISION LIQUID HANDLING



An excellent Combination

The new C Line Syringe from
HAMILTON dedicated to your
CTC  SYSTEMS



HAMILTON

The New Generation of HAMILTON Syringes for CTC ANALYTICS

Together with Hamilton Bonaduz (Switzerland) CTC Analytics has developed a new microliter syringe. The so called C- Line is the next generation of the CTC Mikroliter syringe. It has the following main features.


- Zero carry over
- New Teflon/polymer mixture as sealing material.
Better temperature stability
Surface is less porous.
- Flange with colour code indicating the needle gauge and internal diameter
gauge 22 blue

Syringes for GC PAL (C- Line)

Hamilton P/N	Hamilton Syringe description	CTC P/N
203185	7701.2 CTC	SyrC L1.2-26P-AS
203189	75 N CTC (26S/AS)	SyrC L5-26S-AS
203072	701 N CTC (26S/2)	SyrC L10-26S-2
203205	701 N CTC (26S/AS)	SyrC L10-26S-AS
203361	701 N CTC (23S/AS)	SyrC L10-23-AS
203362	701 N CTC (23S/26S/AS)	SyrC L10-23S/26S-AS
203043	1702 N CTC (26S/AS)	SyrC G25-26S-AS
203076	1710 N CTC (26S/AS)	SyrC G100-26S-AS
203078	1725 N CTC (26/AS)	SyrC G250-26-AS

Replacement Plunger for GC PAL (C- Line)

Hamilton P/N	Hamilton Syringe description	CTC P/N
203245	1702 Plunger CTC 10pcs/pack	Plg G25
203246	1710 Plunger CTC 10pcs/pack	Plg G100
203249	1725 Plunger CTC 10pcs/pack	Plg G250



Postfach 1207 · 52374 Langerwehe
Telefon 024 23 / 40493-0 · Fax 024 23 / 40493-49
E-Mail: info@cs-chromatographie.de