

HAMILTON Cross Reference List HPLC Columns



HAMILTON	PRP-1 (L21)	PRP-3 (L21)	PRP-X100	PRP-X110	PRP-X110S	PRP-X200	PRP-X300	PRP-X400	PRP-X500 (L23)
Waters	Shodex RSpak DS-613		IC-Pak Anion			Shodex RSpak DC-613	IC-Pak Ion Exclusion		Protein-Pak Q 8HR
Dionex	IonPac NS1 IonPac NS1-5um		IonPac AS7		IonPac AS4A IonPac AS4A-SC IonPac AS12A IonPac AS14	IonPac CS10 IonPac CS11	IonPac ICE-AS1		ProPac PA1 ProPac SAX-10
Vydac	259VHP								300VHP
Polymer Labs	PLRP-S 100 A	PLRP-S 300 A							
Metrohm					Phenom. Star Ion A300 Metrosep A Supp 3 IC	Hyperrez Monovalent IC	Metrosep Org. Acids IC		
Phenomenex	Polymerx RP-1				Star Ion A300		Rezex RFQ Fast Acid		
Bio-Rad							Fast Acid Analysis Col. Aminex Glyphosate Col.		
Varian	Chrompack P300 RP								
MN	Nucleogel RP 100	Nucleogel RP 300	Nucleogel Anion I						Nucleogel Ion 300 OA
Supelco							Supelcogel C-610H		
Transgenomic	RPSep PRX-1				ICSep AN300 IC Sep ANSC ICSep AN2 ICSep AN1-SC		ICSep Ion-300		
Shodex			IC I-524A	IC NI-424	IC IF-424	IC Y-521		Shodex IEC QA-825	
DyChrom	Chemco Gel RP col.								
BioChrom Labs	Hydrocell RP 5T								
Merck	Polyspher RP18		Polyspher IC AN	Polyspher IC AN-:Polyspher IC AN-S		Polyspher OA AC		Polyspher OA HY	
Amersham	Source 5RPC Source 15RPC								
MetaChem	MetaChem pHRP								

HAMILTON	PRP-X600	PRP-X700	PRP-X800	RCX-10	RCX-30	HC-40 Ca	HC-75 H (L17)	HC-75 Ca (L19)	HC-75 Pb (L34)
Waters			IC-Pak Ion Cat. col.				Fast Fruit Juice IC-Pak Ion Exclusion Shodex RSpak DC-613	Sugar-Pak 1 Shodex SC-1011 Shodex SC-1821	Shodex SP0810
Dionex	DNAPac PA-100 ProPac WAX-10		IonPac CS5A	CarbPac PA1	CarbPac PA10				
Vydac	301VHP								
Polymer Labs							PL Hi-Plex H	PL Hi-Plex Ca	PL Hi-Plex Pb
Metrohm	Metrosep Carb 1								
Phenomenex							Rezex RHM Monosac.	Rezex RCM Monosac. Rezex RCU-UPS Sugar Alc.	Rezex RPM Monosac.
Bio-Rad						Aminex HPX-42C	Aminex HPX-87H	Aminex HPX-87C	Aminex HPX-87P
Varian							Chrompack Org. Acids	Chrompack Carbohydrates Ca	Chrompack Carbohydrates Pb
MN							Nucleogel Sugar Ca	Nucleogel Sugar Pb	Nucleogel Sugar Pb
Supelco							Supelcogel H	Supelcogel Ca	Supelcogel Pb
Transgenomic			ICSep CN2				CarboSep CHO-620 CarboSep CHO-820 CarboSep USP L19 Ca CarboSep Coregel-87C	CarboSep CHO-682 Lead CarboSep Coregel-87P	
Shodex	RSpak DC-613		IC YK-421				Sugar SH 1011/1821	Sugar SC1011/1821	Sugar SP0810
Alltech							Brownlee Polypore H	Brownlee Polypore Ca	
Merck			Polyspher IC CA				Polyspher OA KC	Polyspher CH CA	Polyspher CH PB
Amersham									
MetaChem							MetaCarb 87H	MetaCarb 87C	MetaCarb 87P